

Dual Protection
Safety I.V. Catheter

Supercath™ 5



MEDIKIT

Color Code ISO10555-5	Catalog Number	Gauge	Catheter			Flow Rate mL/min
			O.D	I.D	Length	
Orange	SP125-14-31T	14 G	2.1 (mm)	1.6 (mm)	1 1/4 " (31 mm)	272
Medium grey	SP125-16-31T	16 G	1.7 (mm)	1.3 (mm)	1 1/4 " (31 mm)	183
Deep green	SP125-18-31T	18 G	1.3 (mm)	0.9 (mm)	1 1/4 " (31 mm)	100
Pink	SP125-20-31T	20 G	1.1 (mm)	0.7 (mm)	1 1/4 " (31 mm)	55
Deep blue	SP125-22-31T	22 G	0.9 (mm)	0.6 (mm)	1 1/4 " (31 mm)	28
Deep blue	SP125-22-25T	22 G	0.9 (mm)	0.6 (mm)	1 " (25 mm)	31
Yellow	SP125-24-19T	24 G	0.7 (mm)	0.4 (mm)	3/4 " (19 mm)	16
Yellow	SP125-24-14T	24 G	0.7 (mm)	0.4 (mm)	0.56 " (14 mm)	21
Violet	SP125-26-19T	26 G	0.6 (mm)	0.4 (mm)	3/4 " (19 mm)	12
Violet	SP125-26-14T	26 G	0.6 (mm)	0.4 (mm)	0.56 " (14 mm)	14
Medium grey	SP120-16-31T	16 G	1.7 (mm)	1.3 (mm)	1 1/4 " (31 mm)	183
Deep green	SP120-18-50T	18 G	1.3 (mm)	0.9 (mm)	2 " (50 mm)	90
Deep green	SP120-18-31T	18 G	1.3 (mm)	0.9 (mm)	1 1/4 " (31 mm)	100
Pink	SP120-20-50T	20 G	1.1 (mm)	0.7 (mm)	2 " (50 mm)	52
Pink	SP120-20-31T	20 G	1.1 (mm)	0.7 (mm)	1 1/4 " (31 mm)	55
Pink	SP120-20-25T	20 G	1.1 (mm)	0.7 (mm)	1 " (25 mm)	62
Deep blue	SP120-22-31T	22 G	0.9 (mm)	0.6 (mm)	1 1/4 " (31 mm)	28
Deep blue	SP120-22-25T	22 G	0.9 (mm)	0.6 (mm)	1 " (25 mm)	31
Yellow	SP120-24-19T	24 G	0.7 (mm)	0.4 (mm)	3/4 " (19 mm)	16
Orange	SP125-14-31-WT	14 G	2.1 (mm)	1.6 (mm)	1 1/4 " (31 mm)	272
Medium grey	SP125-16-31-WT	16 G	1.7 (mm)	1.3 (mm)	1 1/4 " (31 mm)	183
Deep green	SP125-18-31-WT	18 G	1.3 (mm)	0.9 (mm)	1 1/4 " (31 mm)	100
Pink	SP125-20-31-WT	20 G	1.1 (mm)	0.7 (mm)	1 1/4 " (31 mm)	55
Deep blue	SP125-22-31-WT	22 G	0.9 (mm)	0.6 (mm)	1 1/4 " (31 mm)	28
Deep blue	SP125-22-25-WT	22 G	0.9 (mm)	0.6 (mm)	1 " (25 mm)	31
Yellow	SP125-24-19-WT	24 G	0.7 (mm)	0.4 (mm)	3/4 " (19 mm)	16

Yellow	SP125-24-14-WT	24 G	0.7 (mm)	0.4 (mm)	0.56 " (14 mm)	21
Violet	SP125-26-19-WT	26 G	0.6 (mm)	0.4 (mm)	3/4 " (19 mm)	12
Violet	SP125-26-14-WT	26 G	0.6 (mm)	0.4 (mm)	0.56 " (14 mm)	14
Medium grey	SP120-16-31-WT	16 G	1.7 (mm)	1.3 (mm)	1 1/4 " (31 mm)	183
Deep green	SP120-18-50-WT	18 G	1.3 (mm)	0.9 (mm)	2 " (50 mm)	90
Deep green	SP120-18-31-WT	18 G	1.3 (mm)	0.9 (mm)	1 1/4 " (31 mm)	100
Pink	SP120-20-50-WT	20 G	1.1 (mm)	0.7 (mm)	2 " (50 mm)	52
Pink	SP120-20-31-WT	20 G	1.1 (mm)	0.7 (mm)	1 1/4 " (31 mm)	55
Pink	SP120-20-25-WT	20 G	1.1 (mm)	0.7 (mm)	1 " (25 mm)	62
Deep blue	SP120-22-31-WT	22 G	0.9 (mm)	0.6 (mm)	1 1/4 " (31 mm)	28
Deep blue	SP120-22-25-WT	22 G	0.9 (mm)	0.6 (mm)	1 " (25 mm)	31
Yellow	SP120-24-19-WT	24 G	0.7 (mm)	0.4 (mm)	3/4 " (19 mm)	16

50pcs/box

※Please refer to the package insert for complete instruction for use.

※Please contact us from our website for other specifications and availability.



Medikit Co.,Ltd.

1-13-2, Yushima, Bunkyo-ku, Tokyo 113-0034, Japan

TEL:+81-3-3839-0492 | FAX:+81-3-3839-0644 | www.medikit.co.jp/english/

Manufacturer :

Togo Medikit Co., Ltd.

17148-6, Aza Kamekawa, Oaza Hichiya,

Hyuga City, Miyazaki Prefecture 883-0062, Japan

Authorized Representative in the European Community :

Emergo Europe B.V.

Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands

CE 0197

2023.2
20220214RY83-2